

Drugs of Abuse

Cross Reactivity Guide

CEDIA® Benzodiazepine HS 200 ng/mL cutoff

Positive Compounds	Trade Name	Concentration (ng/mL)
7-NH2-Flunitrazepam		210
7-NH2-Nitrazepam		250
a-OH-Alprazolam		120
a-OH-Triazolam		130
Alprazolam	Tafil, Valeans, Xanax, Xandor	100
Alprazolam Glucuronide		200
Bromazepam	Lectopam	200
Chlordiazepoxide	Libritabs, Librium, Sereen	1,300
Clobazam		340
Clonazepam	Klonopin	290
Clorazepate	Gen-XENE, Tranxene	270
Delorazepam		110
Demoxepam		1,100
Desalkylflurazepam		120
Diazepam	Valium, Valrelease, Vazepam	130
Estazolam	ProSom	90
Flunitrazepam	Rohypnol	190
Flurazepam	Dalmane, Durapam	110
Halazepam	Paxipam	200
Lorazepam	Ativan	180
Lorazepam Glucuronide		450
Lormetazepam		150
Medazepam	Nobrium	170
NH2-Clonazepam		210
Nimetazepam		37.5
Nitrazepam		200
Nordiazepam		120
Oxazepam	Serax	160
Oxazepam Glucuronide		800
Prazepam	Centrax	180
Temazepam	Restoril, Razepam, Temaz	220
Temazepam Glucuronide		800
Triazolam	Halcion	100

Negative Compounds	Trade Name	Concentration (ng/mL)
Oxaprozin	Daypro	10,000